



PERFECT FOR EVERYTHING PERFECTA

Perfecta can do it all: from rejuvenation by correcting red and brown skin discoloration, like Aesthetica, to eliminating virtually all vascular and pigmented lesions, like the Platinum. The Perfecta features advanced micro-pulse technology, multiple handpieces for treating various spot sizes, and our patented Dynamic Cooling Device™ (DCD™) cooling system, which all but eliminates patient downtime by preventing purpura, or bruising. And like all Candela products, Perfecta is backed by Candela's unmatched training, marketing and service—including a financing program that's right for your practice.

One laser for all these treatment conditions:

- Skin Rejuvenation
- Diffuse Redness
- Pigmentation
- Freckles
- Sun Spots
- Rosacea
- Acne
- Facial Veins
- Leg Veins
- Pigmented Lesions
- Wrinkles
- Poikiloderma of Civatte
- Angiomas
- Hemangiomas
- Port Wine Stains
- Scars and Stretch Marks
- Warts
- Venous Lakes
- Psoriasis
- Benign Gynecological Vascular Lesions
- Other Benign Cutaneous Vascular Lesions



PERFECT FOR EVERYTHING

See for yourself Perfecta's unmatched versatility.

PRE AND POST-TREATMENT RESULTS



Brown Pigmentation




Courtesy of Taro Kono, M.D.




Brown Pigmentation



Courtesy of Taro Kono, M.D.



Periorbital Wrinkles



Courtesy of Stephen W. Eubanks, M.D.



Redness



Courtesy of David A. Laub, M.D.



Leg Veins



Courtesy of David B. Vasily M.D.



Port Wine Stain




Courtesy of Gerald N. Goldberg, M.D.




Hemangioma



Courtesy of Ross S. Kaplan, M.D.



Nose Veins



Courtesy of Pablo Bioxeda, M.D.



Poikiloderma of Civatte



Courtesy of Stephen W. Eubanks, M.D.



Acne Redness



Courtesy of Wichai Hongcharu, M.D.




Rosacea



Courtesy of J. Charlie Finn, M.D.



Rosacea



Courtesy of Ross S. Kaplan, M.D.



Port Wine Stain



Courtesy of Fine Lines Cosmetic Center



Venous Lake



Courtesy of Fine Lines Cosmetic Center



Scar



Courtesy of Jeffrey S. Dover, M.D.




Warts



Courtesy of Catherine Smadja, M.D.

Redesigned for Maximum Ease of Use

Perfecta has been designed for the ultimate ease of use. Each console is sleeker and smaller than ever before for ease of maneuverability from treatment room to treatment room. Yet, it features a large 10" screen with touch-screen technology for a user-friendly interface. Onscreen graphics are easy to read and provide clear, simple-to-follow treatment guidelines.



Complies with FDA performance standards for laser products except for deviations pursuant to laser Notice No. 50 dated July 26, 2002



Vbeam PERFECTA SPECIFICATIONS

Laser Type	Pulsed Dye
Wavelength	595 nm
Pulse	0.45-40 ms
Spot Sizes (diameter) with maximum energy	3 mm - 40 J/cm ² 5 mm - 30 J/cm ² 7 mm - 20 J/cm ² 10 mm - 10 J/cm ² 12 mm - 7 J/cm ² 3x10 mm - 25 J/cm ² 7 mm PL - 15 J/cm ² 10 mm PL - 10 J/cm ²
Beam Delivery	Lens-coupled 3 m optical fiber with handpiece
Pulse Control	Fingerswitch, Footswitch
Dimensions	43" H x 16" W x 30" D** (109 cm x 40 cm x 76 cm)
Weight	290 lbs (132 kg)
Electrical	230 V 50/60 Hz, single phase, 18A

**39" (99 cm) D with dye cartridge

DYNAMIC COOLING DEVICE

Integrated controls, cryogen container and handpiece with distance gauge

Cryogen	HFC 134a
DCD Spray	User adjustable range:
Duration	20 - 100 ms
DCD Delay	User adjustable range:
Duration	10 - 100 ms

Specifications subject to change without notice.

Vbeam is a registered trademark of Candela Corporation. Dynamic Cooling Device and DCD are trademarks.
©2005 Candela Corporation. All rights reserved. Printed in the USA. 11/05 0920-20-0056 Rev. 04



CANDELA

Candela Corporation
530 Boston Post Road
Wayland, MA 01778 USA
Phone: (508) 358-7637
Fax: (508) 358-5569
Toll Free: (800) 821-2013
www.candelalaser.com