




# Youlaser MT

## *Vaginal Accessories*

**Diva Tight**  
*Regain your intimacy*

## Diva Tight for Vaginal Treatments

### Introduction





Diva Tight is an accessory to be used in combination with YouLaser MT for the treatment of vaginal mucosa's by GaAs 1540 nm and CO2 10600 nm lasers. It uses the Mixed Technology of YouLaser MT to perform an exclusive mixed treatment of the vaginal mucosa by non-ablative and ablative wavelengths in the same time. Diva Tight uses a new scanner specifically developed for the treatment the vaginal walls. The scanner has a wide range of different probes for treating the vagina inside and for the external vulvar side. The new full-angle probe emits the laser energy in a 360 degree angle in one time.

The graduated scale, together with the dedicated speculum, allows a complete control on the treatment uniformity. Faster and less invasive procedures are now possible even for the more severe atrophies.

The scanner can be easily adapted to all patient's needs such as vaginal atrophy, urinary incontinence, vulvar treatment and is also available as upgrade kit for all YouLaser MT already installed.



## Diva Tight Kits

	Code	Description
		<b>Vaginal kit composed by:</b>
		N.1 Scanner (disinfectable)
		N.1 Vaginal probe with pyramidal reflector (not sterile, sterilizable)
	OAM001855.00	N.1 Lens with mounting for vaginal probe
		N.1 Mechanical coupler
		<b>In addition, one of the following:</b>
	MBM006195.00	N.1 Speculum, treatment length 9.5 cm (not sterile, sterilizable)
		<b>Or</b>
	MBM005918.01	N.1 Speculum, treatment length 13 cm (not sterile, sterilizable)
	OAM001734.00	<b>Vulvar Kit</b> Composed by: N.1 Lens with mounting for vulvar tip N.1 Vulvar tip (not sterile, sterilizable)

*Limited warranty: One year for all scanners, optics excluded. Optics and serializable parts are not covered by warranty*



## Diva Tight Spare Parts

	Code	Description
	MBM006198.00	Pyramidal reflector (not sterile, sterilizable)
	OAM001730.00	Lens with mounting for vaginal probe
	OBM000453	Lens for vaginal probe
	MAM000571.00	Vaginal probe with pyramidal reflector (not sterile, sterilizable)
	MBM006195.00	Speculum 9.5 cm treatment length (not sterile, sterilizable)
	MBM005918.01	Speculum 13 cm treatment length (not sterile, sterilizable)
	MAM000572.00	Vulvar tip (not sterile, sterilizable)
	OAM001731.00	Lens with mounting for vulvar tip
	OBM000452	Lens for vulvar tip
	OAM001732.00	Diva Tight Scanner with finger-switch (disinfectable)
	OAM001733.00	Diva Tight Scanner, no finger-switch (disinfectable) for all Youlaser MT with s/n up to xxxxxxx-1115

Limited warranty: One year for all scanners, optics excluded. Optics and serializable parts are not covered by warranty



For more information

[www.quantasystem.com](http://www.quantasystem.com)

[sales@quantasystem.com](mailto:sales@quantasystem.com)

**Diva Tight**  
*Regain your intimacy*